

# Evaluation of ASR for Automatic EEG Artifact Removal

**Chi-Yuan Chang, Sheng-Hsiou Hsu, Luca Pion-Tonachini, and Tzyy-Ping Jung**


Swartz Center for Computational Neuroscience,  
Institute for Neural Computation,  
University of California San Diego

Chang et al., *IEEE Transactions on Biomedical Engineering*, 67(4): 1114 - 21, 2019.


# Artifact Subspace Reconstruction (ASR)




Christian A Kothe




1. ASR selects and concatenate “clean” sections of data from raw data and use them as the reference data.



2. ASR uses the **reference data** and a **user-defined cutoff parameter  $k$**  to determine threshold for artifact recognition.



3. Apply ASR and the threshold back on raw data and reconstruct artifact contaminated portion.



C. Chang, et al., 2019  
C. Chang, et al., 2018  
C.A. Kothe, et al., 2016 (U.S. patent)  
T Mullen, et al., 2015



# Real-time EEG Source-mapping Toolbox (REST)



Luca Pion-Tonachini



Live streaming multi-channel data (LSL)

## Preprocessing:

- Re-reference
- Band-pass filtering (IIR)
- Data cleaning (ASR)

## Source separation:


- Online Recursive ICA

## Source-level analysis:

- Source localization
- Source classification

## Real-time applications:

- Automatic artifact removal
- Brain-state monitoring




Pion-Tonachini & Hsu, *IEEE EMBC*, 2015


Hsu et al., *IEEE EMBC*, 2015

Hsu et al., *IEEE TBME*, 2016

Pion-Tonachini & Hsu et al., *IEEE EMBC*, 2018

Download: <https://github.com/goodshawn12/REST>





- Real-time EEG Source-mapping Toolbox (REST):  
<https://github.com/goodshawn12/REST>


L Pion-Tonachini, et al., 2018

L Pion-Tonachini, et al., 2015




# Challenges: lack of validation

1. How to select ASR threshold for removing artifact components?



2. How effective is the ASR on removing artifacts? Are brain signals affected by ASR?


Artifact-contaminated data




ASR




Artifact-corrected data



# Dataset: EEG during Simulating Driving in a Real Car




# Data preprocessing and analysis



- Remove bad channel(s) in preprocessing
  - ASR will include bad channels when concatenating clean sections.
- Apply ASR with different values for the cutoff parameter  $k$
- Apply ICA to the ASR-corrected data

# Results: EEG Signals Modified by ASR



Chang et al., *IEEE Transactions on Biomedical Engineering*, 67(4): 1114 - 21, 2019.







# Changes in IC Scalp Maps

1. Apply ICA to EEG data with and without ASR cleaning.
2. Obtain ICs' scalp maps
  - Projection of IC activities on the scalp.
3. Calculate the correlation coefficient between scalp maps with and without ASR cleaning.

Non-ASR Cleaning





Chang et al., *IEEE Transactions on Biomedical Engineering*, 67(4): 1114 - 21, 2019.

# Changes in IC Scalp Maps



## Cross-subject Result

1. Use *ICLabel* to classify ICs into **Brain, Muscle, Eye, and Other**.
2. Calculate the dipolarity of each IC.
3. Count the number of retained ICs ( $r > 0.8$ ) in each class after ASR cleaning.




- Classifier: *ICLabel* (Luca Pion-Tonachini, 2019)  
<https://labeling.ucsd.edu/tutorial/overview>

# Changes in IC Activities in a Sample Subject




# Changes in IC activities (a Group Study)

(A) Power of source activities



(B) The percentage of retained power





# An Experiment Paradigm for an Online Evaluation

- Cued-artifact experiment: Motion-related artifacts

Artifact-free Oddball Experiment

Eye Closed  
Resting

180 sec



Audio Oddball  
Experiment

180 sec

**Target artifacts:**  
**Turning head**  
**Nodding**  
**Walking**

Block 1 & 2

Eye  
Open  
Resting



Head  
Turning  
@1 Hz



Eye  
Open  
Resting



Nodding  
@1 Hz



Eye  
Open  
Resting



Eye Open  
Standing




Walking  
@1.5 Hz

80 sec

**Audio oddball experiment is implemented in each session.**




# Event-Related Potential (ERP) in an Auditory Oddball Experiment




# Recover P300 at Cz with ASR during Head Nodding

After FIR band pass filter [1 50]



After ASR with different cutoff parameter



# Summary

- ASR is an **automatic**, **online-capable**, **artifact-removal** method.
- ASR with a mild threshold ( $k = 100$ ) could effectively remove large-amplitude artifacts.
- Previously suggested threshold ( $k = 5-7$ ) could remove up to 50% of brain signal's power.
- This study suggested that a **sensible ASR cutoff parameter  $k$**  might be between **20 to 30**.





